

beta Actin protein (STJ505545)
STJ505545**GENERAL INFORMATION**

Product Type	Proteins
Short Description	beta Actin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60
Gene Symbol	ACTB
Uniprot ID	ACTB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 325-375 on human Actin, cytoplasmic 1 protein.
Immunogen Region	325-375
Specificity	
Immunogen Sequence	