

Deformed Epidermal Autoregulatory Factor 1 protein (STJ505547)

STJ505547

GENERAL INFORMATION

Product Type	Proteins
Short Description	Deformed Epidermal Autoregulatory Factor 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54006
Gene Symbol	Deaf1
Uniprot ID	DEAF1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 500-565 on mouse Deformed epidermal autoregulatory factor 1 homolog protein.
Immunogen Region	500-565
Specificity	
Immunogen	
Sequence	