

Deformed Epidermal Autoregulatory Factor 1 protein (STJ505548)

STJ505548

GENERAL INFORMATION

Product Type	Proteins
Short Description	Deformed Epidermal Autoregulatory Factor 1-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54006
Gene Symbol	Deaf1
Uniprot ID	DEAF1_MOUSE
Immunogen	Purified protein.
Immunogen Region	
Specificity	
Immunogen Sequence	