

THAP1 protein (STJ505549)

STJ505549

GENERAL INFORMATION

Product Type	Proteins
Short Description	THAP1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55145
Gene Symbol	THAP1
Uniprot ID	THAP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 150-210 on human THAP domain-containing protein 1.
Immunogen Region	150-210
Specificity	
Immunogen Sequence	