

TSG101 protein (STJ505553)

STJ505553

GENERAL INFORMATION

Product Type	Proteins
Short Description	TSG101-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7251
Gene Symbol	TSG101
Uniprot ID	TS101_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 330-390 on Tumor susceptibility gene 101 protein.
Immunogen Region	330-390
Specificity	
Immunogen Sequence	