

CD300 protein (STJ505555)
STJ505555**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CD300-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146722
Gene Symbol	CD300LF
Uniprot ID	CLM1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 230-290 on human CMRF35-like molecule 1/CD300LF.
Immunogen Region	230-290
Specificity	
Immunogen Sequence	