

CD300 protein (STJ505556)
STJ505556**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CD300-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84152
Gene Symbol	PPP1R1B
Uniprot ID	PPP1B_HUMAN
Immunogen	Purified protein.
Immunogen Region	
Specificity	
Immunogen Sequence	