

CARD15 protein (STJ505561)

STJ505561

GENERAL INFORMATION

Product Type	Proteins
Short Description	CARD15-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Nod2
Uniprot ID	NOD2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 700-780 on mouse Nucleotide-binding oligomerization domain-containing protein 2.
Immunogen Region	700-780
Specificity	
Immunogen Sequence	