

EGR2 protein (STJ505563)

STJ505563

GENERAL INFORMATION

Product Type	Proteins
Short Description	EGR2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13654
Gene Symbol	Egr2
Uniprot ID	EGR2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 60-150 on mouse E3 SUMO-protein ligase EGR2.
Immunogen Region	60-150
Specificity	
Immunogen Sequence	