

FGD3 protein (STJ505565)

STJ505565

GENERAL INFORMATION

Product Type	Proteins
Short Description	FGD3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89846
Gene Symbol	FGD3
Uniprot ID	FGD3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-320 on human FYVE, RhoGEF and PH domain-containing protein 3.
Immunogen Region	250-320
Specificity	
Immunogen Sequence	