

Ferritin protein (STJ505569)

STJ505569

GENERAL INFORMATION

Product Type	Proteins
Short Description	Ferritin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000 IHC: 1:50-1:150 ICC: 1:50-1:150
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2495
Gene Symbol	FTH1
Uniprot ID	FRIH_HUMAN
Immunogen	Full length Ferritin Protein.
Immunogen Region	
Specificity	
Immunogen Sequence	