

Cell Surface Glycoprotein MUC18 protein (STJ505575)

STJ505575

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Proteins |
| Short Description | Cell Surface Glycoprotein MUC18-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 78967 |
| Gene Symbol | Mcam |
| Uniprot ID | MUC18_RAT |
| Immunogen | Synthetic peptide taken within amino acid region 330-353 on rat Cell Surface Glycoprotein MUC18 protein. |
| Immunogen Region | 330-353 |
| Specificity | |
| Immunogen Sequence | |