

Aromatase protein (STJ505577)

STJ505577

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aromatase-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1588
Gene Symbol	CYP19A1
Uniprot ID	CP19A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Aromatase protein.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	