

Cleavage stimulation factor subunit 1 protein (STJ505579)

STJ505579

GENERAL INFORMATION

Product Type	Proteins
Short Description	Cleavage stimulation factor subunit 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1477
Gene Symbol	CSTF1
Uniprot ID	CSTF1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 200-250 on human Cleavage stimulation factor subunit 1 protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	