

HIPK1 protein (STJ505581)

STJ505581

GENERAL INFORMATION

Product Type	Proteins
Short Description	HIPK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	204851
Gene Symbol	HIPK1
Uniprot ID	HIPK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-750 on human Homeodomain-interacting protein kinase 1 protein.
Immunogen Region	700-750
Specificity	
Immunogen Sequence	