

Indoleamine 2 3-dioxygenase 1 protein (STJ505583)

STJ505583

GENERAL INFORMATION

Product Type	Proteins
Short Description	Indoleamine 2 3-dioxygenase 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3620
Gene Symbol	IDO1
Uniprot ID	I23O1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-150 on human Indoleamine 2, 3-dioxygenase 1 protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	