

Indoleamine 2 3-dioxygenase 2 protein (STJ505586)

STJ505586

GENERAL INFORMATION

Product Type	Proteins
Short Description	Indoleamine 2 3-dioxygenase 2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	169355
Gene Symbol	IDO2
Uniprot ID	I23O2_HUMAN
Immunogen	Purified Protein
Immunogen Region	
Specificity	
Immunogen Sequence	