

Occludin protein (STJ505587)

STJ505587

GENERAL INFORMATION

Product Type	Proteins
Short Description	Occludin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83497
Gene Symbol	Ocln
Uniprot ID	OCLN_RAT
Immunogen	Synthetic peptide taken within amino acid region 450-523 on rat Occludin protein.
Immunogen Region	450-523
Specificity	
Immunogen Sequence	