

HIPK3 protein (STJ505591)

STJ505591

GENERAL INFORMATION

Product Type	Proteins
Short Description	HIPK3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10114
Gene Symbol	HIPK3
Uniprot ID	HIPK3_HUMAN
Immunogen	Synthetic peptide taken within amino aci region 860-910 on human Homeodomain-interacting protein kinase 3 protein.
Immunogen Region	860-910
Specificity	
Immunogen Sequence	