

## Loquacious protein (STJ505598)

STJ505598

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | Loquacious-protein if for use in research applications. |
| <b>Applications</b>      | WB  |
| <b>Host/Source</b>       |   |
| <b>Reactivity</b>        |   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      | WB: 1:2, 500-1:5, 000   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                  |
|---------------------------|------------------|
| <b>Gene ID</b>            |                  |
| <b>Gene Symbol</b>        |                  |
| <b>Uniprot ID</b>         |                  |
| <b>Immunogen</b>          | Purified Protein |
| <b>Immunogen Region</b>   |                  |
| <b>Specificity</b>        |                  |
| <b>Immunogen Sequence</b> |                  |