

AMPK alpha 2 protein (STJ505599)

STJ505599

GENERAL INFORMATION

Product Type	Proteins
Short Description	AMPK alpha 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	78975
Gene Symbol	Prkaa2
Uniprot ID	AAPK2_RAT
Immunogen	Synthetic peptide taken within amino acid region 480-530 on rat 5 NA-AMP-activated protein kinase catalytic subunit alpha-2 protein.
Immunogen Region	480-530
Specificity	
Immunogen Sequence	