

GRK6 protein (STJ505602)

STJ505602

GENERAL INFORMATION

Product Type	Proteins
Short Description	GRK6-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2870
Gene Symbol	GRK6
Uniprot ID	GRK6_HUMAN
Immunogen	Purified Protein
Immunogen Region	
Specificity	
Immunogen Sequence	