

TNIK protein (STJ505603)

STJ505603

GENERAL INFORMATION

Product Type	Proteins
Short Description	TNIK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23043
Gene Symbol	TNIK
Uniprot ID	TNIK_HUMAN
Immunogen	Sythetic peptide taken within amino acid region 1250-1300 on human TRAF2 and NCK-interacting protein kinase protein.
Immunogen Region	1250-1300
Specificity	
Immunogen Sequence	