

CD68 protein (STJ505605)

STJ505605

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD68-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12514
Gene Symbol	Cd68
Uniprot ID	CD68_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 240-290 on mouse CD68/Macrosialin protein.
Immunogen Region	240-290
Specificity	
Immunogen Sequence	