

RBM4 protein (Phospho) (STJ505607)

STJ505607

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-RBM4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5936
Gene Symbol	RBM4
Uniprot ID	RBM4_HUMAN
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 300-350 on human RNA binding motif protein 4.
Immunogen Region	300-350
Specificity	
Immunogen Sequence	