

Angiotensin-like 4 protein (STJ505612)

STJ505612

GENERAL INFORMATION

Product Type	Proteins
Short Description	Angiotensin-like 4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57875
Gene Symbol	Angptl4
Uniprot ID	ANGL4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse Angiotensin-related protein 4.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	