

Aurora A protein (STJ505614)

STJ505614

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aurora A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 350-397 on human Aurora Kinase A protein.
Immunogen Region	350-397
Specificity	
Immunogen Sequence	