

Aurora B protein (STJ505616)

STJ505616

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aurora B-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9212
Gene Symbol	AURKB
Uniprot ID	AURKB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-344 on human Aurora Kinase B protein.
Immunogen Region	300-344
Specificity	
Immunogen Sequence	