

VGLUT1 protein (STJ505618)
STJ505618**GENERAL INFORMATION**

Product Type	Proteins
Short Description	VGLUT1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	72961
Gene Symbol	Slc17a7
Uniprot ID	VGLU1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1-50 on mouse Vesicular Glutamate transporter 1 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	