

PACAP protein (STJ505620)

STJ505620

GENERAL INFORMATION

Product Type	Proteins
Short Description	PACAP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116
Gene Symbol	ADCYAP1
Uniprot ID	PACA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 125-176 on human Pituitary adenylate cyclase-activating polypeptide protein.
Immunogen Region	125-176
Specificity	
Immunogen	
Sequence	