

14-3-3 sigma protein (STJ505622)

STJ505622

GENERAL INFORMATION

Product Type	Proteins
Short Description	14-3-3 sigma-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2810
Gene Symbol	SFN
Uniprot ID	1433S_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 150-220 on human 14-3-3 protein sigma/STRATIFIN/SFN protein.
Immunogen Region	150-220
Specificity	
Immunogen Sequence	