

LRRK1 protein (STJ505624)

STJ505624

GENERAL INFORMATION

Product Type	Proteins
Short Description	LRRK1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79705
Gene Symbol	LRRK1
Uniprot ID	LRRK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1950-2015 on human Leucine-rich repeat serine/threonine-protein kinase 1 protein.
Immunogen Region	1950-2015
Specificity	
Immunogen	
Sequence	