

VPS36 protein (STJ505628)
STJ505628**GENERAL INFORMATION**

Product Type	Proteins
Short Description	VPS36-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51028
Gene Symbol	VPS36
Uniprot ID	VPS36_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-310 on human Vacuolar Protein-Sorting-associated protein 36 protein.
Immunogen Region	250-310
Specificity	
Immunogen	
Sequence	