

Lamin B1 antibody-Nuclear Envelope Marker protein (STJ505630)

STJ505630

GENERAL INFORMATION

Product Type	Proteins
Short Description	Lamin B1 antibody-Nuclear Envelope Marker-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4001
Gene Symbol	LMNB1
Uniprot ID	LMNB1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-500 on human Lamin-B1 protein.
Immunogen Region	450-500
Specificity	
Immunogen Sequence	