

Lamin B1 antibody-Nuclear Envelope Marker protein (STJ505631)

STJ505631

GENERAL INFORMATION

Product Type	Proteins
Short Description	Lamin B1 antibody-Nuclear Envelope Marker-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4001
Gene Symbol	LMNB1
Uniprot ID	LMNB1_HUMAN
Immunogen	Purified protein.
Immunogen Region	
Specificity	
Immunogen Sequence	