

## IRF8 protein (STJ505632)

STJ505632

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	IRF8-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	15900
<b>Gene Symbol</b>	Irf8
<b>Uniprot ID</b>	IRF8_MOUSE
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 350-400 on mouse Interferon Regulatory Factor 8 protein.
<b>Immunogen Region</b>	350-400
<b>Specificity</b>	
<b>Immunogen Sequence</b>	