

IRF8 protein (STJ505633)
STJ505633**GENERAL INFORMATION**

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | IRF8-protein if for use in research applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | WB: 1:2, 500-1:5, 000 |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------|
| Gene ID | 15900 |
| Gene Symbol | Irf8 |
| Uniprot ID | IRF8_MOUSE |
| Immunogen | Purified protein. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |