

CC2D1A protein (STJ505635)
STJ505635**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CC2D1A-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54862
Gene Symbol	CC2D1A
Uniprot ID	C2D1A_HUMAN
Immunogen	Purified protein.
Immunogen Region	
Specificity	
Immunogen Sequence	