

CD134/OX40L receptor protein (STJ505636)

STJ505636

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD134/OX40L receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22163
Gene Symbol	Tnfrsf4
Uniprot ID	TNR4_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 250-300 on mouse Tumor Necrosis Factor Receptor Super Family member 4 protein.
Immunogen Region	250-300
Specificity	
Immunogen Sequence	