

MeCP2 isoform 1 protein (STJ505640)

STJ505640

GENERAL INFORMATION

Product Type	Proteins
Short Description	MeCP2 isoform 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	17257
Gene Symbol	Mecp2
Uniprot ID	MECP2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1-50 on mouse Methyl-CpG-binding Protein 2 isoform 1 protein.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	