

STAT6 protein (Phospho) (STJ505642)

STJ505642

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-STAT6-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6778
Gene Symbol	STAT6
Uniprot ID	STAT6_HUMAN
Immunogen	Purified protein.
Immunogen Region	
Specificity	
Immunogen Sequence	