

KMT6/EZH2 protein (STJ505646)
STJ505646**GENERAL INFORMATION**

Product Type	Proteins
Short Description	KMT6/EZH2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified protein.
Immunogen Region	
Specificity	
Immunogen Sequence	