

DNER protein (STJ505647)

STJ505647

GENERAL INFORMATION

Product Type	Proteins
Short Description	DNER-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	92737
Gene Symbol	DNER
Uniprot ID	DNER_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-737 on human Delta and Notch-like Epidermal growth factor-related receptorprotein.
Immunogen Region	700-737
Specificity	
Immunogen	
Sequence	