

FLT3/CD135 protein (STJ505649)

STJ505649

GENERAL INFORMATION

Product Type	Proteins
Short Description	FLT3/CD135-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2322
Gene Symbol	FLT3
Uniprot ID	FLT3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-950 on human Receptor-type tyrosine-protein kinase FLT3 protein.
Immunogen Region	900-950
Specificity	
Immunogen Sequence	