

TYK2 protein (STJ505653)

STJ505653

GENERAL INFORMATION

Product Type	Proteins
Short Description	TYK2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1150-1186 on rat Non-receptor tyrosine-protein kinase TYK2 protein.
Immunogen Region	1150-1186
Specificity	
Immunogen Sequence	