

Twist protein (STJ505655)

STJ505655

GENERAL INFORMATION

Product Type	Proteins
Short Description	Twist-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7291
Gene Symbol	TWIST1
Uniprot ID	TWIST1_HUMAN
Immunogen	Synthetic peptide take within amino acid region 50-100 on human Twist-related protein 1.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	