

WEE2 protein (STJ505663)
STJ505663**GENERAL INFORMATION**

Product Type	Proteins
Short Description	WEE2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	381759
Gene Symbol	Wee2
Uniprot ID	WEE2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 60-110 on mouse Wee1-like protein kinase isoform 2.
Immunogen Region	60-110
Specificity	
Immunogen Sequence	