

Brk protein (STJ505665)
STJ505665**GENERAL INFORMATION**

Product Type	Proteins
Short Description	Brk-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	20459
Gene Symbol	Ptk6
Uniprot ID	PTK6_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 400-450 on mouse Protein-Tyrosine Kinase 6 isoform.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	