

IRAKM protein (STJ505669)

STJ505669

GENERAL INFORMATION

Product Type	Proteins
Short Description	IRAKM-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11213
Gene Symbol	IRAK3
Uniprot ID	IRAK3_HUMAN
Immunogen	Interleukin-1 receptor-associated kinase 3
Immunogen Region	
Specificity	
Immunogen Sequence	