

## IRAKM protein (STJ505670)

STJ505670

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	IRAKM-protein if for use in research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: 1:2, 500-1:5, 000
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	11213
<b>Gene Symbol</b>	IRAK3
<b>Uniprot ID</b>	IRAK3_HUMAN
<b>Immunogen</b>	Purified Interleukin-1 receptor-associated kinase 3 protein.
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	